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DESCRIZIONE

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RIFERIMENTI

- ¹Strand V. et al., A multicenter, randomized controlled trial comparing a single intra-articular injection of Gel-200, a new cross-linked formulation of hyaluronic acid, to phosphate buffered saline for treatment of osteoarthritis of the knee. Osteoarthritis Cartilage. 20: 350-6, 2012.
- ²Strand V. et al., Effectiveness and safety of a multicenter extension and retreatment trial of Gel-200 in patients with knee osteoarthritis. Cartilage. 3: 297-304, 2012.

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HyLink® INSTRUCTIONS FOR USE

DESCRIPTION

- HyLink® is a sterile, transparent and viscoelastic hydrogel composed of cross-linked hyaluronate, a derivative of highly purified sodium hyaluronate (hyaluronan) extracted from bacterial fermentation.
- Hyaluronan is a polysaccharide containing repeating disaccharide units of D-glucuronic acid and N-acetyl-D-glucosamine. In HyLink®, strands of hyaluronan are bound to each other via dimers of cinnamic acid resulting in increased viscoelasticity.
- HyLink® is delivered in a single-use, prefilled disposable glass syringe.
- This prefilled syringe is composed of a rubber piston [butyl rubber: latex free], rubber tip cap [butyl rubber: latex free], finger grip and plunger rod and is packaged in a molded plastic A-PET film blister with a Tyvek® lid.
- Each prefilled syringe with 3 mL of HyLink® contains:

Cross-linked Hyaluronate	30.0 mg
Sodium Chloride	24.3 mg
Disodium Phosphate Dodecahydrate	0.89 mg
Sodium Dihydrogen Phosphate Dihydrate	1.93 mg
Water for Injection	q.s. to 3.0 mL

INDICATIONS

- Idiopathic osteoarthritis (OA) of the knee.

INTENDED PURPOSE

- HyLink® is indicated for the treatment of pain in idiopathic OA of the knee and is intended to improve the pain by supplementing the viscoelastic hydrogel in OA of the knee. The duration of the effect of an injection is 26 weeks.

INTENDED PATIENT POPULATION

- Patients with idiopathic OA of the knee.

CONTRAINDICATIONS

- Do not administer HyLink® to patients with known hypersensitivity (allergy) to HyLink® or sodium hyaluronate preparations.
- Do not administer HyLink® to patients with severe allergic or autoimmune diseases.

DOSAGE

- HyLink® is a single dose preparation and should only be injected once per treatment course. The dose is 3 mL (one syringe) per knee joint intra-articularly.

ADMINISTRATION

- HyLink® should only be injected by an authorized physician or by an authorized person in accordance with local legislation.
- No additional training is required.
- Check the use-by date indicated on the package.
- Do not use after the use-by date.
- HyLink® must be injected at room temperature. HyLink® must be removed from the refrigerator, where applicable, and kept at room temperature for approximately 30-60 minutes prior to injection.
- Strict aseptic administration technique must be followed.
- Remove joint effusion, if present, through an 18-20 G needle before injecting HyLink®. Maintain needle placement in the joint while disconnecting the syringe used to relieve joint effusion. Discard the syringe containing the removed joint effusion. The same syringe should not be used for both removing effusion and injecting HyLink®.
- Injection of subcutaneous lidocaine or similar local anesthetic may be performed prior to injection of HyLink®.
- Peel off the blister Tyvek® lid from the blister package and remove the syringe.
- Carefully remove the tip cap of the syringe and aseptically attach the syringe to an 18-20 G needle. To ensure a tight seal and to prevent leakage during administration (or injection), secure the needle tightly while firmly holding the luer lock. If effusion was previously removed, connect the syringe to the needle already placed in the joint. Twist the tip cap before pulling it off to minimize product leakage.
- Inject HyLink® into the knee joint through the needle using aseptic injection technique.
- Inject the full, 3 mL of HyLink®, into the knee. If treatment is being administered to both knees, use a separate syringe of HyLink® for each knee.
- For knee joint cavity injection, insert the needle horizontally or slightly downwards in the joint in the gap between the patella and femur. A slight resistance may often be felt as the needle goes through the joint capsule.
- Discard the syringe, the needle and any unused HyLink® in accordance with the institution's policy for infectious waste disposal.

WARNINGS

- HyLink® should be used with caution in patients with infections or skin diseases in the area of the injection site.
- HyLink® should be used with caution if patients have hypersensitivity to any medications.
- HyLink® should be used with caution in patients with systemic bleeding disorders.
- HyLink® should be used with caution in patients with lymphatic or venous stasis in the leg to be injected.
- Do not concomitantly use disinfectants containing quaternary ammonium salts (e.g. benzalkonium chloride) for skin preparation because sodium hyaluronate can precipitate in their presence. Disinfectants such as povidone-iodine are recommended.
- Do not inject HyLink® intravascularly.

PRECAUTIONS

- Strict aseptic administration technique must be followed.
- Remove joint effusion, if present, before injecting HyLink®.
- The safety and performance of the use of HyLink® in joints other than the knee and for conditions other than OA have not been established.
- The safety and performance of the use of HyLink® concomitantly with other intra-articular injectable have not been established.
- Use caution when injecting HyLink® into patients who are allergic to cinnamons.
- The safety and performance of the use of HyLink® in severely inflamed knee joints have not been established.
- Do not inject HyLink® extra-articularly or into the synovial tissue and capsule.
- STERILE CONTENTS: The prefilled syringe is intended for single use. The contents of the syringe must be used immediately once the packaging is opened. Discard any unused HyLink®.
- Do not use HyLink® if the blister package has been opened or damaged, or if there are cracks or breakage in the prefilled syringe.
- Severe joint inflammation in the osteoarthritic knee should be alleviated before injecting HyLink®, as the inflammation may be aggravated by injection.
- Do not use HyLink® after the use-by date indicated on the package.
- As with any invasive joint procedure, it is recommended that the patient avoid any strenuous activities (such as jogging, tennis, other active sports, heavy lifting) and prolonged weight-bearing activities (such as standing for more than one hour) within 48 hours after treatment with the intra-articular injection of HyLink®.

USE IN SPECIFIC POPULATIONS

- Pregnancy: The safety and performance of HyLink® have not been established in pregnant women.
- Lactation: It is not known if HyLink® is excreted in human milk. The

safety and performance of HyLink® have not been established in lactating women.

- Pediatrics: The safety and performance of HyLink® have not been demonstrated in pediatric patients.
- OA severity: The safety and performance of HyLink® have not been demonstrated in patients with a Kellgren-Lawrence Score of 4.
- Body Weight: The safety and performance of HyLink® have not been demonstrated in obese patients (Body Mass Index >35).
- OA Type: Performance of HyLink® has not been established in post-traumatic OA patients.

PERFORMANCE

Improvement in pain and function of the knee joint over 3 to 13 weeks and in pain over 26 weeks after the first injection.

- Based on the results in clinical investigations of SI-6606/01¹, HyLink® demonstrated performance by week 3 after the first injection and it lasted through week 13. WOMAC^(a) VAS^(b) pain subscores were used for the primary outcome measure and it was significant at both week 13 and over 3 to 13 weeks after the first injection (*P* <0.05). HyLink® showed statistically significant improvements from the baseline compared to that in PBS^(c) groups over 3 to 13 weeks after the first injection in WOMAC VAS function subscores, WOMAC VAS total score and physician global evaluation (*P* <0.05). The rate of OMERACT-OARSI strict responders^(d) in the HyLink® group was statistically significant over the PBS group from 6 to 13 weeks after the first injection.
- Based on the re-analyses results of Gel/1133 clinical investigation, a single injection of HyLink® demonstrated significantly more effective than PBS for the reduction of pain in the idiopathic OA of the knee at 26 weeks after the first injection both in WOMAC VAS pain subscores and 50-foot walk test VAS Pain score (*P* = 0.0061).

^(a)WOMAC: Western Ontario and McMaster Universities Osteoarthritis Index.

^(b)VAS: visual analog scale.

^(c)PBS: phosphate buffered saline.

^(d)OMERACT-OARSI strict responders: defined by Outcome Measures in Rheumatoid Arthritis Clinical Trials-Osteoarthritis Research Society International as improvements from baseline in WOMAC pain or physical function subscores ≥50% with absolute changes ≥20 mm (termed "strict responders").

CLINICAL BENEFITS

Improvement in pain and function of the knee joint: over 3 to 13 weeks and in pain over 26 weeks after the first injection.

- Based on the SI-6606/01¹ clinical investigation, HyLink® showed statistically significant improvements from the baseline in HyLink® compared to that in PBS groups over 3 to 13 weeks after the first injection in WOMAC pain subscores, WOMAC function subscores, WOMAC total score and physician global evaluation (*P* <0.05).
- Based on the re-analyses results of Gel/1133 clinical investigation, HyLink® resulted in a statistically significant improvement from the baseline WOMAC pain subscores compared with PBS at 26 weeks after the first injection (*P* = 0.0061).

UNDESIRABLE SIDE-EFFECTS

Reported Undesirable Side-effects:

The following undesirable side-effects have been reported during clinical investigations.

Investigation No.	Undesirable Side-Effects
SI-6606/01 ¹ (N=249)	Joint swelling (14.1%), Joint effusion (11.2%), Arthralgia (7.6%), Injection site pain (2.0%), Joint stiffness (0.8%), Muscular weakness (0.8%), Dizziness (0.8%), Erythema (0.8%), Effusion (0.4%), Injection site bruising (0.4%), Injection site erythema (0.4%), Swelling (0.4%), Alanine aminotransferase increased (0.4%), White blood cell count increased (0.4%), Back pain (0.4%), Muscle spasms (0.4%), Synovitis (0.4%), Tension headache (0.4%), Rash (0.4%), Rash pruritic (0.4%), Hypertension (0.4%)
Gel/1132 ² (N=199)	Arthralgia (9.0%), Joint swelling (8.5%), Joint effusion (7.0%), Injection site pain (1.5%), Blood alkaline phosphatase increased (0.5%), Blood urea increased (0.5%), Burning sensation (0.5%), Contusion (0.5%), Effusion (0.5%), Injection site bruising (0.5%), Injection site reaction (0.5%), Joint warmth (0.5%), Migraine (0.5%), Muscular weakness (0.5%), Swelling (0.5%), White blood cell count decreased (0.5%)
Gel/1133 Double-blind Phase (N=404)	Arthralgia (3.2%), Joint swelling (1.7%), Joint effusion (1.0%), Injection site joint pain (1.0%), Injection site joint swelling (0.5%), Injection site effusion (0.2%), Joint crepitation (0.2%), Joint stiffness (0.2%), Oedema peripheral (0.2%), Urticaria (0.2%)
Gel/1133 Open-label Phase (N=644)	Arthralgia (3.1%), Injection site joint pain (1.4%), Headache (0.5%), Injection site joint swelling (0.5%), Joint swelling (0.5%), Presyncope (0.3%), Procedural dizziness (0.3%), Anxiety (0.2%), Dizziness (0.2%), Joint effusion (0.2%), Joint range of motion decreased (0.2%)

There were neither serious adverse events nor pseudoseptic reactions related to HyLink® injections in the clinical investigations. Transient pain, swelling, and/or effusion of the treated knee joint may occur after intra-articular injection of HyLink®. These events are usually resolved on their own or with conservative treatment. Injection site infections associated with the injection procedure were reported during the post-market experience.

Potential Undesirable Side-effects:

According to post-market experience of other sodium hyaluronate preparations, anaphylactic/anaphylactoid reactions accompanied by transient hypotension (sudden drop in blood pressure), have been rarely reported worldwide and all resolved either spontaneously or after conservative treatment.

Note: If any serious incident occurred in relation to HyLink®, it should be reported to the Authorized Representative, Market Distributor, Manufacturer or Competent authority.

HOW SUPPLIED

- HyLink® is supplied in a 3 mL, disposable, prefilled glass syringe containing 3 mL of HyLink®.
- The contents of the syringe are sterile.
- The product is not made with natural rubber latex.

USE-BY DATE AND STORAGE

- The use-by date is indicated on the labelling and the package.
- Store in the original package at a temperature of 1°C – 25°C.
- DO NOT FREEZE.

REFERENCES

¹ Strand V. et al., A multicenter, randomized controlled trial comparing a single intra-articular injection of Gel-200, a new cross-linked formulation of hyaluronic acid, to phosphate buffered saline for treatment of osteoarthritis of the knee. Osteoarthritis Cartilage. 20: 350-6, 2012.

² Strand V. et al., Effectiveness and safety of a multicenter extension and retreatment trial of Gel-200 in patients with knee osteoarthritis. Cartilage. 3: 297-304, 2012.

Emergo Europe B.V.
Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

Importer :

TRB CHEMEDICA S.R.L.
Via Mauro Macchi 27, 20124 Milano, Italy
email: contact@trbchemedica.it

SEIKAGAKU CORPORATION

6-1, Marunouchi 1-chome, Chiyoda-ku
Tokyo 100-0005, Japan
Telephone: +81-3-5220-8974
www.seikagaku.co.jp/en

© Trademark of SEIKAGAKU CORPORATION, Japan

0123

Summary of safety and clinical performance

Summary of safety and clinical performance (SSCP) for HyLink is available in the European database on medical devices (Eudamed). The Eudamed public website: <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 4987541HyLinkW5 is used to search and find the SSCP for HyLink in the Eudamed.

LIST OF SYMBOLS

Symbols	Descriptions
	Authorized representative in the European Community
	Importer
	Manufacturer
	Unique Device Identifier
	Sterilized using steam or dry heat. Sterile fluid path.
	Single sterile barrier system with protective packaging outside
	Do not re-use
	Consult instructions for use
	Caution
	Do not use if the package is opened or damaged
	Temperature limit
	Medical Device
	Batch code
	Use-by date
	Month and year of manufacture
	CE mark (compliance with Regulation (EU) 2017/745)
	Reference